

# A multiplexed star electrode array for session-independent silent speech recognition

Sudharshan Sundaramahalingam Supervisor: Dr Arsam Shiraz



# Silent speech recognition

Speech recognition without the use of audio signals, through other correlated bio-signals. Allows for detection of silent/mouthed speech.

Is as intuitive and fast as speech but is also:



Private and convenient in public.



Unaffected by environmental noise



Accessible to people with difficulty vocalizing

High impact applications such as:



**Defense** ~ 500B USD market by 2030



**Human-Al interaction** ~ 100B USD market by 2030



# **EMG-based silent speech recognition**

Electromyography (EMG) signals from facial muscles is a common modality for silent speech recognition.

### **Typical Architecture:**

Electrodes

Signal Acquisition

Digital Filters


Word Classifier



State of the art
AlterEgo (2019), 91% accuracy
over a 30-word vocabulary

### Replicated state-of-the-art with a custom electrode brace





System performance requires precise electrode positions, accuracy reduces to 40% after electrode shift.

### Research Question:


Investigating the impact of electrode alignment with facial muscles to maximise signal-tonoise ratio (SNR) in silent speech interfaces.



# **Design and Development**

Develop an 8-channel electrode array capable of actively switching the orientation of its electrodes.

### System Design



### System Development



**Main Control Board** 



**Multiplexer Board** 





Full Setup

(top) multiplexed, (bottom) control

Slide 4 of 7




# **Experimental Procedure and Development**

Identify electrode directions which maximize SNR, and investigate variability in directions due to electrode shift.

### **Methodology**

- 1. Record one minute of noise in a relaxed state.
- Record a phonetically balanced passage that ensures a variety of speech muscle activations are sampled.
- 3. Calculate SNR using mean-squared value of noise and recorded signal.
- 4. Repeat for each electrode direction (0, 45, 90, 135) + all directions + no directions selected + control array.
- 5. Repeat across 2 sessions to sample variability in the presence of electrode shift.



**Detailed view of electrode array:** showing electrode directions and channels



**User wearing electrodes:** Arrows indicate orientation compared to figure above.

Electrodes adhered with bio-safe double-sided tape.




# **Key Results**

Optimally oriented electrodes consistently outperform average control electrode SNR of 10 dB.

The optimal orientation of the electrode **changes significantly** between sessions.

Multiplexing can offer significant SNR improvements.





### **Conclusion**

- A novel approach to optimizing SNR for silent speech recognition through electrode alignment is investigated.
- A multiplexing system enabling configurable electrode orientation is developed.
- Results show evidence that dynamic electrode alignment has the potential to boost SNR significantly compared to a static electrode array.
- 4. System can be improved with automatic electrode alignment for user-administrable silent speech interfaces that are robust to electrode placement.

